

Procedure following Discharge of Defibrillator

If possible request a replacement set of pads from paramedics

A

If pads obtained:

Replace new pads tray in defibrillator and run through check procedure

If unit functioning and rescue ready replace in cabinet (if not follow B instructions)

Contact Clerk for replacement rescue kit

B

If pads not obtained:

Put Sign in defibrillator cabinet informing out of use and location of nearest alternative

Put defibrillator in a safe place (NOT back in cabinet)

Contact Clerk and inform of defibrillator discharge

Clerk will inform EMAS to de-activate CPAD until informed site is rescue ready

(If Clerk unavailable 0115 9193417 or 07966806507)

Clerk will order new pads

(then follow A instructions)